



Dental Burs

Ultrasound Scalers

Radicular Posts

Teeth Whitening

# CATALOGUE

**2022-23**





# Welcome

We are pleased to share with our dear clients our **Axis Dental Catalogue 2022-23**, in which we present our wide range of products.

Since the foundation of our company in Switzerland in 1998, we have worked towards the goal of continuous innovation that leads us to offer the highest quality and efficiency in all our products. This goal of finding excellence is the key behind the high satisfaction of dentists, laboratory technicians and, of course, dental patients.

Our wish to continuously improve has led us in recent years to a process of digital transformation of our company. After a process of rebranding, we are launching a brand new WEBSITEsite and our very first online store, an window of new opportunities to offer our products worldwide in the most agile way.

We also continue to serve our clients through our team of sales representatives located in Switzerland, Spain, Portugal and Bulgaria. We are well know by our clientele from our dedicated service and attention to detail that allows to observe the clients´ needs.

In this catalogue you will find, in addition to the references of our more than **1000 modalities of dental burs, ultrasound scalers, radicular posts and dental whitening treatments**, useful information for the dentist on product maintenance and conservation of our dental instruments. Because the quality and durability of our products is what defines us and on which the trust and loyalty of our customers is based upon.

Thanks for choosing us.

Jean-Paul Berclaz      Dominique Rahm

Partners  
[Axis Dental](#)

**AXIS Dental Sarl**  
Ch. du Cloalet 4  
1023 Crissier - Swizerland  
Tel: +41 21 633 0 633  
Email: [mail@axis-dental.ch](mailto:mail@axis-dental.ch)

[www.axis-dental.com](http://www.axis-dental.com)

Design and layout: [www.laurasalesa.com](http://www.laurasalesa.com)

Our *Guarantee* lies  
upon our experience.

## CONTENTS



Our company



Guarantee



International network



Products

## DENTAL BURS



### Introduction and features


#### Clinic

Diamond burs  
Tungsten Carbide Burs  
Ceramics  
Polishers

#### Laboratory

Diamond burs  
Tungsten Carbide Burs  
Polishers

12



## RADICULAR POSTS



### Introduction and features

Fibra de vidrio

192



## TEETH WHITENING



### Introduction and features

Clinic  
At home

216





## INTRODUCTION



**Axis Dental is a company specialized in the manufacture, selling and distribution of dental burs, ultrasound scalers, radicular posts and teeth whitening solutions.**

### Phylosophy

Our philosophy is based on the manufacture of precise and efficient instruments of high quality with high durability for the long term benefit of the dentist

#### **Personalized treatment: dedicated and professional**

Our business goal is to win the trust of dental professionals through the performance of a team of professionals that have been trained to offer a personalised service to clients.




**Our clients appreciate the high quality and durability of Axis Dental products.**

**Dominique Rahm**  
*Sales and Marketing, Axis Dental*



### Mission



Axis Dental seeks excellence in the quality of its products and services. The professionalism of the Axis Dental's staff meets the needs and expectations of our clients.



*Our company*





## INTRODUCTION



**Created in 1998 in Switzerland, Axis Dental aims to offer professionals in the dental sector the best solutions that facilitate their work thanks to materials of the latest generation.**

### Quality and experience

Our dental instruments are developed and produced in our manufacture based in Switzerland. This is the warranty of the high quality of all our products.

Using the latest technology and equipment, we monitor each phase of the production process to ensure the best results. All Axis Dental products meet the most strict European standards.



We consider it is essential to offer our product honesty. Our trust on our products' quality makes obvious for the client that we have a good product.



### Commitment

Axis Dental seeks excellence in the quality of its products and service. Precision, quality and service are the values that make us evolve.

Axis Dental has earned the trust of professionals thanks to its personalized service, its human factor and the high-quality of its products.

Our guarantee is based on our goal for continuous improvement and adaptation to current processes and regulations.

### International quality management

Our medical devices are produced in accordance with international standards ISO 13485.

### European Community Certificate

All our medical products have CE recognition and comply with the provisions of Directive 93 /42 /EEC.



## Guarantee





## INTRODUCTION



**Our company has its own manufacturing laboratory in Switzerland, as well as delegations in Spain, Portugal and Bulgaria. We are also present through distributors in various countries.**

### International distribution

Axis Dental has an international sales network capable of providing quality support to our customers and partners.

At Axis Dental there is always an expert near you. We have a large team of sales technicians, so that our products are accessible to all clinics and laboratories in our distribution areas.

The professionalism of the Axis Dental staff can competently meet the needs and expectations of our clients..

*International network*

### Close to you


Our excellence in customer relationship relationship is based on 4 pillars:

- ✓ We are customer service orientated
- ✓ We reward your loyalty
- ✓ We customize Kits according to your needs
- ✓ We offer training on our products

### Personalized kits

Our solutions include a wide range of references in dental burs for both clinical and laboratory work. All of them are made with the best materials, to provide them with efficiency, precision, cutting speed that allow them to create perfect geometries with zero vibrations.

We can create a personalized kit that fits the needs and work dynamics of the dentist, both at the clinic as well as at the laboratory.



### Professional training

As part of our commitment to the development of current and future professionals, Axis Dental in Spain has developed AD Training, a training program adapted to advances and new work techniques in the field of dentistry.

Our workshops, specialized in the areas of dental restoration and aesthetics, have a theoretical-practical component that allows attendees to reinforce the concepts taught by renowned professionals from around the world.

Discover our training offer available to develop your professional skills at the highest level.








# DENTAL BURS

AXIS  
dental  
SWISS QUALITY




*An essential  
tool in  
Dentistry*





# DENTAL BURS



## An essential tool in Dentistry

Dental burs occupy an important place in the basic instruments of Dentistry. Its correct choice and adaptation to the rotary instrument will determine its long-term efficiency.

Very often these instruments are chosen empirically. The study of the different substances and materials on which we work and the results of different clinical studies, allow us to choose the burs that best adapt.

### Normas de uso

- Enough cooling is needed through water spray application in order to prevent deterioration of the tooth and its surroundings.
- The spray must be distributed over the length of the instrument.
- A deflection of the jet could cause injuries of thermal origin.
- Additional cooling is necessary for instruments with a length greater than 19 mm or with a diameter of the active part greater than 2 mm.
- Immediately remove damaged instruments
- Take note of the prescribed working pressure of 0.3 - 2N (30 - 200p) for FG instruments and 2 - 5N (200 - 500p) for HP instruments, work without pressure.
- The turbines, contra-angles and handpieces must be in perfect condition for the optimum performance of the instrument.
- Insert instruments carefully and without forcing them in the

training mode. The instruments must be inserted into the forceps until the end of their cylindrical part.

- During work, avoid tilting the instrument and turning it on the edge.
- You must respect the recommended rotation speeds indicated on the packaging.
- Bring the instrument to its recommended rotation speed outside the patient's mouth (for the prosthetist bring the rotation to the maximum speed before applying the instrument to the object to be worked on).
- Any worn or damaged instruments must be disposed of in accordance with the waste disposal regulations in force in the country of use.
- Any serious incident that occurs in relation to the instrument must be reported to the manufacturer and/or to the competent authority of the Member State in which the user is established.

## >>> Dental burs family



At Clinic

### DIAMOND

- Standard
- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental



At Laboratory

### AT LABORATORY

- Diamond
- Tungsten Carbide Burs
- Polishers
- Diverse



### TUNGSTEN CARBIDE BURS

- Standard
- Excavation
- Crowns and bridges
- Finishing 12 sheets
- Finishing 30 sheets
- Root canal
- Surgery




### CERAMIC

- Ceramic for soft tissue



### POLISHERS

- Brushes
- Universal
- Composites
- Ceramics
- Stones



### Applications



Maxillofacial Surgery



Prophylaxis



Implants



Removal of old restorations



Orthodontics



Crown separation



Crown preparation



Root planing



Cavity preparation



Treatment of fillings



Root canal preparation





The main characteristics of the dental burs are:

- The strict choice of the material with which they are made.
- The precision and rigor of manufacturing.
- The effectiveness and comfort of the product when working.
- The lifetime of the instrument in the medium-long term.






# DENTAL BURS



On our **WEBSITE** you will find all our bur references with detailed information to be able to choose the products that suit the needs of the professional.


**VISIT THROUGH THE QR.**



**VISIT OUR**

# WEBSITESITE

**How to place an order:**



## 1 With Axis reference code

Adding the diameter to the indicated Axis reference.

**FG F 801.014**

## 2 Through our e-shop

By directly selecting the desired diameter of the tip in the e-shop. The system will automatically suggest the Axis reference.

**FG G 801.027**

## 3 With the ISO number

Adding the diameter of the indicated ISO number.

**806.313.001.524.010**




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

Standard  
R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

**Axis Dental has designed a new generation of instrument, a bur unique with 20% more coverage by Diamond.**



### Introduction to Diamond burs

Diamond instruments are made of a stainless steel handle (monobloc) that is high resistance to deformation. It is perfectly calibrated to have a pure rotation movement.

The diamond is fixed on top of the bur, joint to the bare steel through a special electrodeposition process in three phases:

- A nickel preparation layer that plays a role
- of anti-vibration.
- An anchoring layer, which bonds the Diamond grains together securely.
- A hard layer of Nickel that strongly fixes the grains by Diamond.

They can be manufactured with natural or synthetic Diamond particles. Considering the superiority of natural Diamond particles, Axis Dental only manufactures its products with this option.

The physical properties of natural Diamond are as follows:

- Its high hardness.
- His great purity.
- Its way.
- Its structure.

The active part of the bur is covered with calibrated natural Diamond grains. The great abrasive capacity of this surface, associated with a high rotation speed, allows a very effective abrasion effect on the hard structures of the tooth (enamel and healthy dentin).

However, it is not effective in removing soft tooth structure (carious dentin).



### How to use:

- The bur is fixed on a high speed turbine or contra angle.
- These devices have a water and air cooling system.
- The rotation of the cutter and the Size of the Diamond particles vary according to the desired efficiency required.
- Due to their high efficiency, coarse or super coarse diamond instruments are used to perform any type of dental preparation without increasing the temperature.
- The standard burs are of medium grain particles, allowing a good abrasion and a fine cut.
- Diamond finishing instruments have a fine grain that allows a more precise work while obtaining a perfectly smooth surface.


The power of  
Diamond



# DENTAL BURS




Clinic



DIAMOND TUNGSTEN CARBIDE BURS CERAMICS POLISHERS


## ISO number breakdown explanation

(ISO 6360, 2157, 1797)



$A + B C + D + E + F = \text{ISO-No. } 806.314.111.534.012$

<b>A</b>	Active part's material Example: 806 Diamond / 500 Tungsten Carbide Burs, etc.	<b>E</b>	Active part's finishing Example: 534 bold Diamond grains, 007 helical blades with transverse cuts, etc.
<b>B / C</b>	Handle and total length Example: 314 FG / 204 RA, etc.	<b>F</b>	Nominal size, diameter of the largest active part on scale 1 /10 mm. Example: 012 diameter 0.12 mm, etc.
<b>D</b>	Active part's shape Example: 111 cylider, 001 redondeada, etc.		



## Correspondence of ISO number and Axis Dental's reference

ISO-No. 806.314.111.534.012

FG G 837.012  
Handle Finishing Shape Diameter

ISO	AXIS DENTAL
806	Diamond
314	FG
111	837
534	G
012	012

## Dental burs handles

### Diameter: 1.60 mm

For all types of high-speed turbines or contra-angles (red)

ISO 313 FG S (Friction Grip Short)

- Size of Handle L: 16 mm
- In all cases, in which the working height is limited, a smaller handle allows a more comfortable access.

ISO 314 FG (Friction Grip)


- Size of Handle L: 19 mm
- It is the most common handle.

ISO 315 FG L (Friction Grip Long )

- Size of Handle L: 21 mm

ISO 316 FG XL (Friction Grip Xtra Long )

- Size of Handle L: 25 mm
- This handle's size is optimum for surgery



### Diameter: 2.35 mm

For all types of low speed contra-angles (blue or green)

ISO 204 RA (Right Angle)

- Size of Handle L: 22 mm

ISO 205 RA L (Right Angle Long)

- Size of Handle L: 26 mm


ISO 206 RA XL (Right Angle Xtra Long)

- Size of Handle L: 34 mm

For use in handpieces

ISO 104 HP type (Hand Piece)

- Size of Handle L: 44.5 mm



Fotos proporcionada por BienAir.







# DENTAL BURS

DIAMOND

TUNGSTEN CARBIDE BURS


CERAMICS

POLISHERS




Clinic

AXIS  
dental  
SWISS QUALITY



## Characteristics of Diamond Burs



AXIS  
dental  
SWISS QUALITY

1

### HARDNESS

Natural Diamond is the hardest mineral in the world.

Weight	3.514
Molecular mass	12.01
Hardness	10 MOHS= 8000 H.V

3

### SHAPE

The octahedral shape of natural Diamond particles allows for high Diamond density and regular inter-particle spacing. Both features are necessary for optimal cutting residue removal.

The concentration obtained is between 50 and 60% of the surface, while the optimal theoretical concentration obtained with spheres is 75%.

Diamond particles must be calibrated to ensure an optimal and regular concentration on the instrument. Its equal distribution will guarantee a perfect concentric rotation of the bur.

2

### PURITY

Natural Diamond particles have virtually no magnetic impurities such as nickel, cobalt or iron.


A De Beers study compares the purity of natural Diamond with synthetic.

Comparison Between Natural and Synthetic Diamonds		
	Natural diamond	Synthetic diamond specially treated by electro deposition
Iron impurities Ppm	10-50	1000-7000
Magnetic susceptibility X 10-6 cgs	0-4	110-170

4

### SUPERFICIAL STRUCTURE

The irregular surface of the natural Diamond particles gives it very good retention effects on the nickel bonding layer and improves abrasive behavior thanks to its sharp cutting edges.





# DENTAL BURS



## DIAMOND

TUNGSTEN CARBIDE BURS

CERAMICS

POLISHERS

Calibrated Diamonds particles are mainly distributed in five different sizes.



**Super Thick (SG) - BLACK RING**  
Diamond grain size ~ 150-180µ

- The most abrasive instrument on the market.
- The efficiency is 60% higher than that of a standard cutter.
- They withstand high cutting pressure with a high abrasion capacity.

**Thick (G) - GREEN RING**  
Diamond grain size ~ 125-150 µ

- A relatively abrasive instrument.
- The efficiency is 30% higher than that of a standard milling cutter.
- They withstand high cutting pressure with a high abrasion capacity.

**Standard - NO RING**  
Diamond grain size ~ 106 µ

Good abrasion capacity.

**Thin (F) - RED RING**  
Diamond grain size ~45µ

Smooth surface treated with a good finish.

**Super Thin (C) - YELLOW RING**  
Diamond grain size ~ 15µ

- Smooth surface treated with a good finish.
- Ideal application on composites.

Average life of the abrasive capacity of Diamond instruments based on the materials and substances treated

The average abrasion life of Diamond instruments depends on the work performed and the manner in which they are used. Its wear comes from the detachment of Diamond particles.

The half life of the instruments can decrease with incorrect use.



- Enamel causes normal wear.
- Healthy dentin, cements and composites minimally disturb the instrument.
- The amalgam and gold, rather than wear down the instrument, they clog the space between the Diamond particles.
- Steel, chrome-nickel and chrome-cobalt quickly alter the instrument.
- The tip of a steel post included in the restoration can damage the abrasive surface of the instrument.




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY





## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

### STANDARD

- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

### Standard



AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS



CERAMICS

POLISHERS

### STANDARD

- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

### Standard



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY





## DIAMOND

### STANDARD


- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

### 368




	L mm	3.5	3.5	5	5
806.314.257.524	FG 368	016	018	021	023
806.314.257.504	FG C 368	016	018	021	023
806.314.257.514	FG F 368	016	018	021	023
806.314.257.534	FG G 368	016	018	021	023
806.314.257.544	FG SG 368	016	018	021	023
806.313.257.524	FG S 368	016			023
806.204.257.524	RA 368				023
806.204.257.514	RA F 368				023

### 369





	L mm	5.5
806.314.263.524.023	FG 369	023
806.314.263.504.023	FG C 369	023
806.314.263.514.023	FG F 369	023
806.314.263.534.023	FG G 369	023
806.314.263.544.023	FG SG 369	023



C 379

### 379



	L mm	3	3.5	4.5	4.5
806.314.277.524	FG 379	014	018	021	023
806.314.277.504	FG C 379	014	018	021	023
806.314.277.514	FG F 379	014	018	021	023
806.314.277.534	FG G 379	014	018	021	023
806.314.277.544	FG SG 379		018	021	023




# DENTAL BURS



Clinic


AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental



### 390

Ø 1,6 mm  
19 mm

L mm 3,5

806.314.274.524.016	FG 390	016
806.314.274.504.316	FG C 390	016
806.314.274.514.016	FG F 390	016

### 392

Ø 1,6 mm  
19 mm

L mm 5

806.314.465.524.016	FG 392	016
806.314.465.504.016	FG C 392	016
806.314.465.514.016	FG F 392	016

### 801

Ø 1,6 mm  
19 mm

L mm 0,8 0,9 1 1,2 1,4 1,6 1,8 2,1 2,3 2,7 2,9 3,5

806.314.001.524	FG 801	008 009 010 012 014 016 018 021 023 027 029 035
806.314.001.504	FG C 801	X X X 012 014 X X X 023 X 029 X
806.314.001.514	FG F 801	X 008 010 012 014 016 018 021 023 027 029 035
806.314.001.534	FG G 801	X 009 010 012 014 016 018 021 023 027 029 035
806.314.001.544	FG SG 801	X X X 012 014 016 018 021 023 X X X

Ø 1,6 mm  
19 mm

806.313.001.524	FG S 801	X 009 010 012 014 016 X X X X X
-----------------	----------	---------------------------------

Ø 2,35 mm  
25 mm

806.204.001.524	RA 801	X X X 012 014 016 018 X 023 X 029 X
-----------------	--------	-------------------------------------

### 801L

Ø 1,6 mm  
19 mm

L mm 1 1,2 1,4 1,6 1,8

806.314.697.524	FG 801L	010 012 014 016 018
806.314.697.514	FG F 801L	X X X 016 X
806.314.697.534	FG G 801L	010 012 014 016 018
806.314.697.544	FG SG 801L	X 012 014 016 018






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY






## DIAMOND

### STANDARD




R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

### 802






L mm	2.5	3	3	3	3.5
FG 802	010	012	014	016	018
FG F 802	X	012	X	X	X
FG G 802	010	012	014	016	X

### 805






L mm	1	1	1.2	1.5	1.5	1.5	2
FG 805	009	010	012	014	016	018	023
FG G 805	X	010	012	014	016	018	X

### 806






L mm	2	3	3	3	3
FG 806	010	012	014	016	018
FG G 806	010	012	014	016	018

### 807






L mm	3.5	4	5
FG 807	012	016	018
FG G 807	012	016	018

### 811



L mm	4	6
FG 811	033	040
FG SG 811	033	040

### 813



L mm	1.5	1.6
FG 813	010	014
FG G 813	X	014

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

### STANDARD

- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

## 818






L mm	0,5	1	1
FG 818:	018	033	042
FG G 818:	X	X	042

806.314.042.524

806.314.042.534

## 820



L mm	2
FG 820:	050

806.314.043.524.050

## 822



L mm	2
FG 822:	009

806.314.237.524.009

## 829



L mm	3	3
FG 829:	010	012
FG C 829:	010	012
FG F 829:	010	012

806.314.235.524

806.314.235.504

806.314.235.514

## 829L



L mm	5
FG 829L:	014
FG C 829L:	014
FG F 829L:	014

806.314.236.524.014

806.314.236.504.014

806.314.236.514.014


AXIS  
dental  
SWISS QUALITY



## Clinic



AXIS  
dental  
SWISS QUALITY



## DIAMOND

## TUNGSTEN CARBIDE BURS

CERAMICS

## POLISHERS








# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

### 835

L mm	3	3	4	4	4	5
FG 835	008	009	010	012	014	016
FG C 835	008	X	X	X	X	X
FG F 835	008	X	X	X	X	X
FG G 835	008	009	010	012	014	016
FG SG 835	X	X	010	012	014	X

806.314.107.524	FG 835	008	009	010	012	014	016
806.314.107.504	FG C 835	008	X	X	X	X	X
806.314.107.514	FG F 835	008	X	X	X	X	X
806.314.107.534	FG G 835	008	009	010	012	014	016
806.314.107.544	FG SG 835	X	X	010	012	014	X

806.313.107.524	FG S 835	008	X	010	012	X	X
-----------------	----------	-----	---	-----	-----	---	---

806.204.107.524	RA 835	X	X	X	012	X	X
-----------------	--------	---	---	---	-----	---	---

### 835KR

L mm	3	4	4	4	5
FG 835KR	008	010	012	014	016
FG F 835KR	008	010	012	014	016
FG G 835KR	008	010	012	014	016

806.314.156.524	FG 835KR	008	010	012	014	016
806.314.156.514	FG F 835KR	008	010	012	014	016
806.314.156.534	FG G 835KR	008	010	012	014	016

### 836

L mm	5	6	6	6
FG 836	010	012	014	018
FG C 836	X	012	X	X
FG F 836	010	012	014	X
FG G 836	010	012	014	018
FG SG 836	X	012	014	X

806.314.110.524	FG 836	010	012	014	018
806.314.110.504	FG C 836	X	012	X	X
806.314.110.514	FG F 836	010	012	014	X
806.314.110.534	FG G 836	010	012	014	018
806.314.110.544	FG SG 836	X	012	014	X

806.313.110.524	FG S 836	X	012	X	X
-----------------	----------	---	-----	---	---

### 836KR

L mm	5	6	6	6	6
FG 836KR	010	012	014	016	018
FG F 836KR	X	012	X	016	X
FG G 836KR	010	X	X	X	X
FG SG 836KR	010	X	X	X	X

806.314.157.524	FG 836KR	010	012	014	016	018
806.314.157.514	FG F 836KR	X	012	X	016	X
806.314.157.534	FG G 836KR	010	X	X	X	X
806.314.157.544	FG SG 836KR	010	X	X	X	X




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

### STANDARD

R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

### 837

L mm	7	8	8	8	8
FG 837	010	012	014	016	018
FG F 837	010	012	014	X	X
FG G 837	010	012	014	016	018
FG SG 837	010	012	014	X	X
FG S 837	X	012	X	X	X
RA 837	X	X	014	X	X

### 837KR

L mm	8	8	8
FG 837KR	012	014	016
FG F 837KR	012	014	016
FG G 837KR	012	014	016



G 837

### 837L

L mm	10	10	10	10
FG 837L	012	014	016	018
FG F 837L	012	014	X	X
FG G 837L	012	014	016	018




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

## TUNGSTEN CARBIDE BURS

## CERAMICS

## POLISHERS

### 838



L mm 3 4 4

	806.314.137.524	FG 838	008	010	012
	806.314.137.514	FG F 838	X	X	012
	806.314.137.534	FG G 838	X	010	012
	806.314.137.544	FG SG 838	X	010	012

### 845




L mm 4 4 4

	806.314.170.524	FG 845	012	014	016
	806.314.170.514	FG F 845	012	X	X
	806.314.170.534	FG G 845	012	014	X
	806.314.170.544	FG SG 845	012	X	X

### 845KR



L mm 4 4 4

	806.314.544.524	FG 845KR	018	021	025
	806.314.544.514	FG F 845KR	018	021	025

### 846



L mm 6 6 7

	806.314.171.524	FG 846	012	016	024
	806.314.171.514	FG F 846	X	016	024
	806.314.171.534	FG G 846	012	016	024
	806.314.171.544	FG SG 846	X	016	X



	806.313.171.524	FG S 846	X	016	X
---	-----------------	----------	---	-----	---






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

### 846KR

Technical drawing showing dimensions: Ø 1.6 mm, 19 mm.

	L mm	6	6	6
806.314.545.524	FG 846KR	014	016	018
806.314.545.514	FG F 846KR	X	016	018
806.314.545.534	FG G 846KR	014	016	X
806.314.545.544	FG SG 846KR	014	016	X

### 847

Technical drawing showing dimensions: Ø 1.6 mm, 19 mm.

	L mm	8	8	8	8	8
806.314.172.524	FG 847	012	014	016	018	021
806.314.172.514	FG F 847	012	X	X	018	021
806.314.172.534	FG G 847	012	014	016	018	021
806.314.172.544	FG SG 847	012	014	016	018	021

Technical drawing showing dimensions: Ø 1.6 mm, 16 mm.

806.313.172.524	FG S 847	X	014	X	X	X
-----------------	----------	---	-----	---	---	---

### 847KR

Technical drawing showing dimensions: Ø 1.6 mm, 19 mm.

	L mm	8	8	8	8
806.314.545.524	FG 847KR	014	016	018	023
806.314.545.514	FG F 847KR	X	X	018	023
806.314.545.534	FG G 847KR	X	X	018	023
806.314.545.544	FG SG 847KR	X	X	018	X

### 848

Technical drawing showing dimensions: Ø 1.6 mm, 19 mm.

	L mm	10	10	10	10	10
806.314.173.524	FG 848	012	014	016	018	025
806.314.173.504	FG C 848	X	X	016	X	X
806.314.173.514	FG F 848	012	014	016	X	X
806.314.173.534	FG G 848	012	014	016	018	025
806.314.173.544	FG SG 848	012	014	016	018	X

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD



R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

### 849




Part Number	FG	010	012
806.314.196.524	FG 849	010	012
806.314.196.514	FG F 849	010	012
806.314.196.534	FG G 849	010	012

### 850




Part Number	FG	012	014	016	018	021	023
806.314.199.524	FG 850	012	014	016	018	021	023
806.314.199.504	FG C 850	012	014	016	X	X	X
806.314.199.514	FG F 850	012	014	016	018	021	023
806.314.199.534	FG G 850	012	014	016	018	021	023
806.314.199.544	FG SG 850	012	014	016	018	021	023



Part Number	FG	014	X	X	X	X
806.204.199.524	RA 850	014	X	X	X	X

### 850XL



Part Number	FG	012
806.314.199XL.524.012	FG 850XL	012
806.314.199XL.514.012	FG F 850XL	012
806.314.199XL.534.012	FG G 850XL	012
806.314.199XL.544.012	FG SG 850XL	012

### 851




Part Number	FG	012	016
806.314.219.524	FG 851	012	016
806.314.219.534	FG G 851	012	016



# DENTAL BURS



Clinic



## DIAMOND

- STANDARD
- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

### 852

	L mm	4	6	6
806.314.164.524	FG 852	010	012	015
806.314.164.504	FG C 852	010	012	015
806.314.164.514	FG F 852	010	012	015
806.314.164.534	FG G 852	X	012	015
806.313.164.524	FG S 852	010	X	X

### 855

	L mm	6	6
806.314.197.524	FG 855	012	014
806.314.197.514	FG F 855	012	X
806.314.197.534	FG G 855	012	014

### 856

	L mm	8	8	6	8	8	8
806.314.198.524	FG 856	012	014	015	016	021	023
806.314.198.504	FG C 856	012	014	015	X	X	X
806.314.198.514	FG F 856	012	X	016	018	021	023
806.314.198.534	FG G 856	012	014	015	018	021	023
806.314.198.544	FG SG 856	012	014	015	016	021	023
806.313.198.524	FG S 856	X	014	015	X	X	X






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD



- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

## TUNGSTEN CARBIDE BURS

## CERAMICS

## POLISHERS

### 858





Ø 1,6 mm  
14 mm

	FG 858	012	014	016	018
806.314.165.524					
806.314.165.504					
806.314.165.514					
806.314.165.534					
806.314.165.544					

Ø 1,6 mm  
14 mm

	FG S 858	014	016	018
806.313.165.524				

### 859



Ø 1,6 mm  
10 mm

	FG 859	010	012	014	016	018
806.314.166.524						
806.314.166.504						
806.314.166.514						
806.314.166.534						
806.314.166.544						

Ø 1,6 mm  
16 mm

	FG S 859	010	012	016	018
806.313.166.524					



Ø 2,35 mm  
22 mm

	RA 859	012	016	018
806.204.166.524				



G 858

### 860



Ø 1,6 mm  
10 mm

	FG 860	010	012
806.314.247.524			
806.314.247.504			
806.314.247.514			
806.314.247.534			

Ø 2,35 mm  
22 mm

	RA 860	012
806.204.247.524		

### 861



Ø 1,6 mm  
10 mm

	FG 861	012
806.314.248.524 012		
806.314.248.504 012		
806.314.248.514 012		




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD


- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

## 862



L mm 8 8 8 8

	806.314.249.524	FG 862	010	012	014	016
	806.314.249.504	FG C 862	010	012	014	016
	806.314.249.514	FG F 862	010	012	014	016
	806.314.249.534	FG G 862	010	012	014	016
	806.314.249.544	FG SG 862	X	012	014	016



	806.313.249.524	FG S 862	X	012	X	X
--	-----------------	----------	---	-----	---	---



	806.204.249.524	RA 862	X	012	X	X
	806.204.249.514	RA F 862	X	012	X	X

## 863



L mm 10 10 10

	806.314.250.524	FG 863	012	015	018	
	806.314.250.504	FG C 863	012	015		X
	806.314.250.514	FG F 863	012	015	018	
	806.314.250.534	FG G 863	012	015	018	
	806.314.250.544	FG SG 863	012	015		X



	806.204.250.524	RA 863	X	X	018	
--	-----------------	--------	---	---	-----	--

## 864



L mm 12

	806.314.251.524.014	FG 864	014			
	806.314.251.534.014	FG G 864	014			

## 868K



L mm 8 8

	806.314.222.524	FG 868K	012	016		
	806.314.222.504	FG C 868K	X	016		
	806.314.222.514	FG F 868K	012	016		
	806.314.222.534	FG G 868K	012	016		
	806.314.222.544	FG SG 868K	012	016		

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

## 877



L mm

B

G

806.314.288.524

FG 877

010

012

806.314.288.504

FG C 877

010

X

806.314.288.514

FG F 877

010

012

806.314.288.534

FG G 877

010

012

806.314.288.544

FG SG 877

010

012



806.313.288.524

FG S 877

010

X

## 877K



L mm

G

G

806.314.297.524

FG 877K

014

016

806.314.297.544

FG SG 877K

014

016

## 878



L mm

B

B

B

B

806.314.289.524

FG 878

010

012

014

016

806.314.289.504

FG C 878

010

X

X

X

806.314.289.514

FG F 878

010

012

014

016

806.314.289.534

FG G 878

010

012

014

016

806.314.289.544

FG SG 878

010

012

014

016



806.313.289.524

FG S 878

X

X

014

X

## 878K



L mm

B

B

B

B

B

806.314.298.524

FG 878K

012

014

016

018

021

806.314.298.514

FG F 878K

012

014

016

018

021

806.314.298.534

FG G 878K

012

014

016

018

021

806.314.298.544

FG SG 878K

X

014

016

018

021

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Clinic



## DIAMOND

### STANDARD

R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

## TUNGSTEN CARBIDE BURS

## CERAMICS

## POLISHERS

### 879

Technical drawing showing dimensions:  $\varnothing 1,6 \text{ mm}$  and 10 mm.

Product	FG	012	014	016
806.314.299.524	FG B79	012	014	016
806.314.299.504	FG C 879	012	X	X
806.314.299.514	FG F 879	012	014	016
806.314.299.534	FG G 879	012	014	016
806.314.299.544	FG SG 879	012	014	016

### 879K

Technical drawing showing dimensions:  $\varnothing 1,6 \text{ mm}$  and 10 mm.

Product	FG	012	014	016	018
806.314.299.524	FG B79K	012	014	016	018
806.314.299.514	FG F 879K	012	014	016	018
806.314.299.534	FG G 879K	012	014	016	018
806.314.299.544	FG SG 879K	X	014	016	018

### 880

Technical drawing showing dimensions:  $\varnothing 1,6 \text{ mm}$  and 10 mm.

Product	FG	010	012	014
806.314.140.524	FG 880	010	012	014
806.314.140.514	FG F 880	010	012	014
806.314.140.534	FG G 880	010	012	014
806.314.140.544	FG SG 880	010	X	X

### 881

Technical drawing showing dimensions:  $\varnothing 1,6 \text{ mm}$  and 10 mm.

Product	FG	010	012	014	016	018
806.314.141.524	FG 881	010	012	014	016	018
806.314.141.504	FG C 881	X	X	014	X	X
806.314.141.514	FG F 881	010	012	014	016	018
806.314.141.534	FG G 881	010	012	014	016	X
806.314.141.544	FG SG 881	X	012	014	016	X




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

## TUNGSTEN CARBIDE BURS

## CERAMICS

## POLISHERS

### 882



L mm

10

10

10

806.314.142.524

FG 882

010

012

014

806.314.142.514

FG F 882

010

012

014

806.314.142.534

FG G 882

010

012

014

806.314.142.544

FG SG 882

010

012

014

### 883



L mm

3

806.314.142.524.010

FG 883

010

806.314.142.514.010

FG F 883

010

806.314.142.534.010

FG G 883

010

### 884



L mm

6

806.314.129.524.012

FG 884

012

806.314.129.514.012

FG F 884

012

806.314.129.534.012

FG G 884

012



806.313.129.524.012

FG S 884

012

### 885



L mm

8

806.314.130.524

FG 885

012

806.314.130.514

FG F 885

012

806.314.130.534

FG G 885

012

806.314.130.544

FG SG 885

012

### 886



L mm

10

10

10

806.314.131.524

FG 886

012

014

016

806.314.131.504

FG C 886

012

X

X

806.314.131.514

FG F 886

012

014

016

806.314.131.534

FG G 886

012

014

016

806.314.131.544

FG SG 886

012

X

X

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

## TUNGSTEN CARBIDE BURS

## CERAMICS



## POLISHERS

### 888



Part Number	FG	Color
806.314.496.524.012	FG 888	012
806.314.496.534.012	FG G 888	012

### 889





Part Number	FG	Color
806.314.540.524	FG 889	009 010
806.314.540.504	FG C 889	009 X
806.314.540.514	FG F 889	009 010
806.314.540.534	FG G 889	009 010
806.313.540.524	FG S 889	009 X

### 898



Part Number	FG	Color
806.314.213.524.016	FG 898	016
806.314.213.534.016	FG G 898	016

### 909



Part Number	FG	Color
806.314.068.524	FG 909	033 042 050
806.314.068.534	FG G 909	033 042 X
806.314.068.544	FG SG 909	X 042 X
806.313.068.524	FG S 909	X 042 X
806.204.068.524	RA 909	033 X X






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

- Standard
- R-RD SERIES**
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

### R-RD Series



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS

CERAMICS

POLISHERS

Standard  
R-RD Series  
**MICROPREPARATION**  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

### Micropreparation



830M




830RM




838M



889M





953AM



953M



838M



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS


CERAMICS

POLISHERS

- Standard
- R-RD Series
- MICROPREPARATION**
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental




953M



On our **WEBSITE** you will find all our bur references with detailed information to be able to choose the products that suit the needs of the professional.

**VISIT THROUGH THE QR.**



VISIT OUR

# WEBSITE

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

- Standard
- R-RD Series
- Micropreparation
- VENEERS**
- Scalers guide
- Diamond DLC
- Multilayer
- Periodontology
- Interdental

## VENEERS



834



839



GF 844

### 834

	016	018	021
L mm	6	6	6

806.314.552.524

FG 834

### 839


	012	014	016
L mm	0.2	0.2	0.2

806.314.150.524

FG 839

## Veneers preparation

The guide grooves define the maximum amount of tissue that we remove, thus determining the level of the final preparation. The smallest (016) creates 0.3mm depths on the central incisors and canines, while the longest head size (021) creates 0.5mm depths on the lateral incisors. With the 850 /016 we reduce the buccal structure of the tooth to the predetermined depth.



### GF844

	014	016
L mm	10	10

806.314.150.524

FG GF844






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

Standard  
R-RD Series  
Micropreparation  
Veneers  
**SCALERS GUIDE**  
Diamond DLC  
Multilayer  
Periodontology  
Interdental

With the support of a guide at its end,  
this bur ensures an identical and  
controlled width over the 360 degrees of  
carving the crown.

### Scalers guide



856P



878KP



880P



881P






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY




## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

Standard  
R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
**DIAMOND DLC**  
Multilayer  
Periodontology  
Interdental

**Diamond burs with superficial DLC  
(Diamond-like Carbon) which gives them  
unparalleled durability.**

### Diamond DLC



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY




## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

Standard  
R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
**MULTILAYER**  
Periodontology  
Interdental

Following the latest technological advances, the Multilayer range of Diamond burs is made up of 2 types of grains, which allows working with them on the most resistant material and extends the life of the bur.

### Multilayer



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

- Standard
- R-RD Series
- Micropreparation
- Veneers
- Scalers guide
- Diamond DLC
- Multilayer
- PERIODONTOLOGY**
- Interdental

**Diamond burs for Periodontology have a longer Handle and contra-angle fixation; in this way, the success of the treatments carried out with them is ensured.**

### Periodontology



P368L





P368XL



P868L



P868XL








# DENTAL BURS



Clinic



## DIAMOND

TUNGSTEN CARBIDE BURS CERAMICS POLISHERS

Standard  
R-RD Series  
Micropreparation  
Veneers  
Scalers guide  
Diamond DLC  
Multilayer  
Periodontology  
**INTERDENTAL**

**Diamonds Interdental Burs that require an EVA contra-angle.**

## INTERDENTAL




ORTHO



PROXO


ORTHO	
	5 Åµm / Å" 0.12
	15 Åµm / Å" 0.14
	25 Åµm / Å" 0.16
	45 Åµm / Å" 0.26
	XF ORTHO
	UF ORTHO
	C ORTHO
	F ORTHO

PROXO	
	25 Åµm
	45 Åµm
	60 Åµm
	90 Åµm
	C PROXO
	F PROXO
	PROXO
	G PROXO



**On our WEBSITE you will find all our bur references with detailed information to be able to choose the products that suit the needs of the professional.**

**VISIT THROUGH THE QR.**



**VISIT OUR**

# WEBSITE




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUNGSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

Axis Dental has developed an innovative design for Tungsten Carbide Burs that solves the common problems that existed with this instrument.




1

## CAVITIES REMOVAL

Efficient and comfortable

These burs allow the enamel to be cut without first using a Diamond bur.




2

## REMOVAL OF SOFT RESTORATIONS (amalgams, composites, cements)

Safe and functional

These ranges offer 2 options to suit the way you work.



3

## REMOVAL OF METAL RESTORATIONS



Quality and power

A special carbide alloy that cuts metal very efficiently.

4

## CROWNS AND BRIDGES

Special carbide with crossed blades for better control and comfort over the cutter.



5

## DENTISTRY SECTION

Perfect cutting capacity with optimal control and safety.

Carbide is 10 times more effective than steel.

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUNGSTEN CARBIDE BURS

CERAMICS

POLISHERS

Thanks to the characteristics of the Tungsten Carbide Burs selected by Axis Dental, the cutting effect of our burs is several times higher than a standard steel bur.



**The active tip of these tungsten carbide burs solves 6 problems for the dentist:**

- Avoid shock at the time of contact with matter.
- With our exclusive blade, the cutter enters the material smoothly without resistance and/or impact.
- Eliminates vibrations during work.
- It is a comfortable and pleasant bur to work with where the professional has total control of the bur during the intervention.
- It avoids the rupture of the active part of the strawberr.
- It does not cause unnecessary imbalances in the turbine, preventing costly revisions.

**Tungsten Carbide Burs cutters are the best tool for working on soft materials, for different reasons:**

- They are most effective in these types of situations.
- It manages to leave the surface that the professional needs, smooth or retentive, depending on the treatment.
- They reduce the intervention time of the professional.
- Metals require specific attention in order to avoid heating the area of the intervention and making the patient uncomfortable. Axis Dental's Tungsten Carbide Burs have superior quality and a special design that ensure total success on this treatment.
- Axis Dental always offers you the best choice of each bur for each type of activity.


## Tungsten Carbide Burs

The active part is made of blades whose profile allows an effective cutting effect associated with waste removal.

- The leaves are helical.
- The sheets may have transverse cuts that form small grooves.
- The number of sheets and the space between them may increase or decrease depending on the purpose of the job.
- Few blades separated from each other will produce a good cutting effectiveness and an efficient removal of residues, in exchange for greater fragility and vibration.
- A large number of blades close together will produce a slow and precise finish job without vibration, with less cutting effectiveness and less heat production.

### How to use:

- Mounted on a turbine or high-speed contra-angle for milling hard tooth structures.
- Mounted on a low speed contra-angle for trimming the cavity, expulsion of soft tissues, carious dentin...



## Active extremity ensuring comfort and effectiveness

Axis Dental's round burs are distinguished by an active transverse blade at their extremity, facilitating work without shocks and with the absence of vibrations.

The cutting efficiency of Axis Dental's Tungsten Carbide Burs allows you to work easily on both rough and smooth surfaces and always with the best quality.

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS








POLISHERS

## STANDARD

- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## Standard



1												
												
L mm		0.5	0.6	0.8	1	1.2	1.4	1.6	1.8	2.1	2.3	
												
500.314.001.001		FG 1	005	006	008	010	012	014	016	018	021	023
												
500.313.001.001		FG S 1	X	X	X	010	X	014	X	X	X	X
												
500.316.001.001		FG XL 1	X	X	X	010	X	014	X	018	X	023
												
500.204.001.001		RA 1	005	006	008	010	012	014	016	018	X	023
												
500.104.001.001		HP 1	005	006	008	010	012	014	016	018	X	023

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS






Clinic



## DIAMOND TUGNSTEN CARBIDE BURS CERAMICS POLISHERS

- STANDARD
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

### Standard

<div>2</div> <div></div> <div></div> <div></div> <div></div> <div></div>	<div></div>							
	L mm	0.5	0.9	1.1	1.2	1.4	1.6	1.7
 500.314.010.001	FG 2	006	008	010	012	014	016	018




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

## STANDARD

- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## Standard


### 21



500.314.107.006

FG 21

010 012 014 016



500.204.107.006

RA 21

010 X X X



500.104.107.006

HP 21

010 012 014 X



### 21L



500.314.110.006

FG 21L

010 012



### 21R



500.314.137.006

FG 21R

009 010 012



500.104.137.006

HP 21R

X 010 012



### 23



500.314.168.006

FG 23

010 012 016



500.204.168.006

RA 23

010 012 X



500.104.168.006

HP 23

010 X 016



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

## STANDARD

- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## Standard

### 23L



500.314.171.006

FG 23L 009 012

### 23R



500.314.194.006

FG 23R 010 012 016



500.204.194.006


RA 23R X 012 016



500.104.194.006

HP 23R X 012 016

### 23RL



500.314.196.006

FG 23RL 010 012

### 25R



500.314.238.006

FG 25R 010 012

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS















CERAMICS





POLISHERS







## STANDARD

- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## Standard

<div>31</div> <div></div> <div></div> <div></div>	<div></div>							
	L mm	4.2	4.2	4.2	4.4	4.4		
		500.314.107.007	FG 31	009	010	012	014	016
								
		500.313.107.007	FG S 31	X	X	012	X	X
<div></div> <div></div>								
		500.316.107.007	FG XL 31	X	010	012	014	X
								
	500.204.107.007	RA 31	009	010	012	014	X	
								
	500.104.107.007	HP 31	009	010	012	014	016	

<b>31L</b>					
					
					
					
		500.314.110.007		FG 31L	<b>D10</b> D11

<h1>31R</h1>   <div></div>					
		L mm	3.9	4.2	4.2
 500.314.137.007		FG 31R	005	010	012
					
 500.104.137.007		HP 31R	009	X	012



AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS








CERAMICS





POLISHERS



## STANDARD

- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## Standard

33					
					
		L mm	4.2	4.2	4.4
					
					
					
					
500.314.158.007		FG 33	010	012	016
500.316.168.007		FG XL 33	X	012	016
500.204.168.007		RA 33	010	012	X
500.104.168.007		HP 33	010	012	016

33L					
					
		L mm	5.5	5.5	5.5
					
500.314.171.007		FG 33L	009	010	012
500.104.171.007		HP 33L	X	010	012

33R					
					
		L mm	4.2	4.2	4.4
					
500.314.194.007		FG 33R	010	012	016
500.104.194.007		HP 33R	X	X	016




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- EXCAVATION**
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## EXCAVATION



### 1S

L mm	1	1.2	1.4	1.6	1.8	2.1	2.3
FG 1S	010	012	014	016	018	021	023
RA 1S	010	012	014	016	018	021	023
RA L 1S	010	012	014	016	018	021	023



### Q1

L mm	1	1.2	1.4	1.6	1.8	2.1	2.3
FG Q1	010	012	014	016	018	021	023
RA Q1	010	012	014	016	018	021	023
RA L Q1	010	012	014	016	018	021	023




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUNGSTEN CARBIDE BURS




CERAMICS

POLISHERS

- Standard
- Excavation
- CROWNS AND BRIDGES**
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

Specific milling cutters  
manufactured with a special  
Tungsten Carbide Burs, which  
allows them to deal with high  
efficiency in cutting the metals of  
the Crowns and bridges.

CROWNS AND BRIDGES



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND


TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 SHEETS FINISHING**
- 30 sheets finishing
- Root canal
- Surgery

## 12 SHEETS FINISHING



### 41



500.314.001.072

L mm 1.4 1.8 2.3


FG 41 014 018 023



500.204.277.072

RA 41 014 018 023

### 44E



500.314.277.072

L mm 3.1 3.5 4.2

FG 44E 014 018 023



500.204.277.072

RA 44E X 018 023

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 SHEETS FINISHING**
- 30 sheets finishing
- Root canal
- Surgery

## 12 SHEETS FINISHING

### 46



Ø 1,6 mm  
19 mm

500.314.254.072

FG 46

L mm	3.5	3.5	3.8	4
014				
018				
023				

Ø 2,35 mm  
22 mm

500.204.254.072

RA 46

015			
-----	--	--	--

### 47L



Ø 1,6 mm  
19 mm

500.314.234.072

FG 47L

L mm	4.4	4.4
012		
014		

Ø 2,35 mm  
22 mm

500.204.234.072

RA 47L

012	
-----	--

### 48




Ø 1,6 mm  
19 mm

500.314.246.071

FG 48

L mm	3.6	3.6	3.6
009			
010			
012			

### 48L



Ø 2,35 mm  
22 mm

500.204.249.072.012

RA 48L

L mm	8
012	

### 50



Ø 1,6 mm  
19 mm

500.314.465.041.012

FG 50

L mm	3.2
012	




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 SHEETS FINISHING**
- 30 sheets finishing
- Root canal
- Surgery

## 12 SHEETS FINISHING

### 132



500.314.181.072

L mm 3

FG 132 009

### 133



500.314.217.072

L mm 5.5

FG 133 010

### 135



500.314.217.072

L mm 8

FG 135 014



212RL

### 212




500.314.197.072.012

L mm 5.5

FG 212 012

### 212R



500.314.197.072.010

L mm 6

FG 212R 010

### 212RL



500.314.198.072.016

L mm 9


FG 212RL 016



# DENTAL BURS



Clinic



DIAMOND


TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS


- Standard
- Excavation
- Crowns and bridges
- 12 SHEETS FINISHING**
- 30 sheets finishing
- Root canal
- Surgery

12 SHEETS FINISHING



On our **WEBSITE** you will find all our burs references with detailed information to be able to choose the products that suit the needs of the professional.

VISIT THROUGH THE QR.



VISIT OUR

# WEBSITE




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 SHEETS FINISHING**
- Root canal
- Surgery

## 30 SHEETS FINISHING



U 41



U 44E



U 48



U 50



U 212



U 247








# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 SHEETS FINISHING**
- Root canal
- Surgery

30 SHEETS FINISHING

## U212



500.314.197.102.010

FG-U212 010

L mm 15

## U247



500.314.499.101.012

FG-U247 012

L mm 3.4

**THE MOST AGILE FORMAT TO HAVE EVERYTHING AT HAND.**

**Our solutions include a wide range of references in dental burs for both Clinic and Laboratory work.**

All of them incorporate the best materials providing efficiency, precision, cutting speed, perfect geometries and absence of vibrations.

We can create a kit personalized for you.

PERSONALISED

## DENTAL BURS KITS



AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND


TUGNSTEN CARBIDE BURS

CERAMICS


POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal**
- Surgery

ROOT CANAL



152



500.315.210.295.014

FG L 152 014



500.316.210.295.014

FG XL 152 014

560



500.316.196.336.008

FG XL 560 008




AXIS  
dental  
SWISS QUALITY



# DENTAL BURS



Clinic



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery


## SURGERY



## WE ARE NEAR YOU

Our extensive network of national and international sales representatives is always attentive to the needs of its Clinics and Laboratories in the area.

We are here to help.



## INTERNATIONAL NETWORK




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery

## SURGERY


### 141



500.205.001.251

RA L 141

010 014 018 023 027 031 035 040 050



500.206.001.251

RA XL 141

010 014 018 023 027 031 035 040 050



500.104.001.251

HP 141

010 014 018 023 027 031 035 040 050

### 151



500.315.199.295.014

FG L 151

014



500.316.199.295.014

FG XL 151

014

L mm

11

## SURGERY

### 161



500.315.408.295

FG L 161

012

014

L mm

10

10

### 162A



500.316.408.298.016

FG XL 162A

016

L mm

10



500.205.408.298.016

RA L 162A

016



500.104.408.298.016

HP 162A

016






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND

TUGNSTEN CARBIDE BURS

CERAMICS

POLISHERS

- Standard
- Excavation
- Crowns and bridges
- 12 sheets finishing
- 30 sheets finishing
- Root canal
- Surgery**

## SURGERY

### 162



500.316.408.297.018

FG XL 162 018



500.205.408.297.018

RA L 162 018



500.104.408.297.018

HP 162 018

### 166



500.205.409.297.021

RA L 166 021



500.104.409.297

HP 166 021

## SURGERY

### 167



500.205.410.297

RA L 167 023



500.104.410.297

HP 167 023

### 389



316.494.020.014

FG XL 389 014




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS CERAMIC POLISHERS

**The ideal solution for cutting soft tissue while cauterizing the wound.  
A simple and effective alternative to the electric scalpel.**



## Application

- Clear deep cavities.
- Interradicular granular tissue laceration.
- Papillectomy / remove hyperplastic gingiva.
- Exposure of sections of osseointegrated implants and included teeth.

Before using CSTT, the area should be cleaned and dried using a cotton swab or air syringe. The tip of the CSTT should be moved through the tissue in small brush strokes at a 30-45 degree angle, avoiding forcing the tip of the trimmer into the tissue or cutting too deeply in a single cut.

The CSTT moves in the area in a balanced traversing motion to the required depth. To avoid any loss of tissue, it is important to preserve the epithelium of the marginal gingiva.

Using the CSTT correctly, in all its applications, you will see that the CSTT is not only capable of replacing conventional and electric scalpels, but it can also prevent us from using the retraction cord on many occasions.

## Example



Cutting soft tissue

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS CERAMIC POLISHERS

## CERAMIC



CSTT



CSTT L



## Information to user

The CSTT is used with high turbine rotation (300,000 - 450,000 rpm). When correctly used, the kinetic energy produces coagulation by heat without burning the tissue. When using the CSTT, avoid excessive contact with teeth or bone tissue.

## Cooling


Without water or cooling spray

## Cleaning

Water and detergent. Important: do not use any solvent.

## Sterilization

With dry heat at a maximum temperature of 180°C or with a steam autoclave at 134°C for 18 minutes.






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

**Different solutions of Polishers with additives adapted to any type of material found in the dental field.**



## NEW GENERATION OF POLISHERS

The exceptional polishing power of these instruments adds extra value to the restorations made. The brilliant result that you will achieve with the new generation of Axis Dental Polishers will keep the beauty and perfection of the work carried out in the long term.

This new generation of Polishers makes us feel proud of all the effort made by the company in research and manufacturing. Our success is based on our mission to achieve high-quality products by introducing innovative materials. And our success is rewarded through customer satisfaction.

The professional faces many problems when performing the treatments, some of them unexpected, such as allergies to certain materials. These responsibilities of the healthcare professional are very important to Axis Dental and are part of our commitment to the Dentist.

Replacing the silicone of traditional Polishers with polyurethane manages to solve numerous allergy problems that affect a part of the patients.

## REVOLUTIONS AND PROTOCOL OF USE

We advise respecting the revolutions for the use of the Polishers, thus ensuring maximum efficiency of the same, getting the professional to work faster, the polisher wears less, extends its useful life and therefore saves material for the Clinic.

An excess of revolutions will dull the polisher, the heat will rise immediately, nullifying the effectiveness of the polisher and also the adequate aesthetic result will not be achieved.

Axis Dental professionals make up a team trained to be able to solve all your doubts regarding the appropriate revolutions or any question regarding the correct use of the material. Contact the delegate in your area who will be happy to help you.

## PRODUCT APPLICATION

The polisher is used with a water spray and without exerting pressure so that the dentist will only have to seek contact with the material.



*A new generation of dental polishers*

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Clinic



DIAMOND | TUGSTEN CARBIDE BURS | CERAMIC | POLISHERS

## TYPE AND ABRASION



Universal polishers for amalgam, gold, non-precious metal alloys, titanium.

With polyurethane and silicon carbide

10 Series...	Efecto espejo	Abrasion 60 / 105 microns
20 Series...	Prepolished	Abrasion 20 / 25 microns
30 Series	Mirror effect	Abrasion 2 / 4 microns

### Polishers for composites

With polyurethane and silicon carbide

Serie 40..	Prepolished	Abrasion 45 / 75 microns
Serie 50..	Mirror effect	Abrasion 2 / 4 microns

### Polishers for composites and hybrid Ceramic

With polyurethane and Diamond

Serie X45..	Prepolished	Abrasion 50 / 65 microns
Seria X55..	Mirror brightness	Abrasion 4 / 8 microns

Polishers for porcelain (Ceramic), lithium disilicate, ceramic oxide, porcelain

With polyurethane and Diamond


Serie X70..	Prepolished	Abrasion 40 / 58 microns
Serie X80..	Mirror effect	Abrasion 10 / 20 microns

### Brushes

- PROPHYLAXIS
- Nylon
- Diamond

### Stones

- Green
- Arkansas






# DENTAL BURS



Clinic



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- BRUSHES
- Universal
- Composite
- Ceramic
- Stones

### Brushes / PROPHYLAXIS

## BPH



Ø mm

5

RA BPH

### Brushes / PROPHYLAXIS nylon

## BPNC



Ø mm

6

RA BPNC

### Brushes / PROPHYLAXIS nylon

## BPNN



Ø mm

4

RA BPNN

## BPNP



Ø mm

4

RA BPNP




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUNGSTEN CARBIDE BURS CERAMIC POLISHERS

- BRUSHES
- Universal
- Composite
- Ceramic
- Stones

### Brushes / Diamond

#### BP D C



Ø mm 6

RA BP D C

#### BP D N



Ø mm 4

RA BP D N

#### BP D P



Ø mm 4

RA BP D P


### CHOICE OF SHAPE

#### How to choose the right shape?

When choosing the shape of a polisher, not only must the type of treatment to be carried out be taken into consideration, it is also important to take into account the savings that we can have in the purchase of material.

Within the same family, the life of the product can be much higher. For example, the 4014 shape (disk shape) has a 2x lifespan than a 4006 shape (point shape).

This situation improves even more when we compare the 4012 shape (wheel shape) which has a 2.5 times higher lifespan than the 4006 (point shape).




AXIS  
dental  
SWISS QUALITY



# DENTAL BURS



Clinic



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- UNIVERSAL
- Composite
- Ceramic
- Stones

### UNIVERSAL / Efecto espejo

## 1005



L mm	9.5
Ø mm	6.5




653.204.030.533

RA 1005

## 1006





L mm	7.5
Ø mm	3



653.314.243.533

FG 1006



653.204.243.533

RA 1006

### UNIVERSAL / Efecto espejo

## 1008



L mm	10
Ø mm	5



653.204.300.533

RA 1008

## 1012



L mm	1.5
Ø mm	14



RA 1012




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- UNIVERSAL
- Composite
- Ceramic
- Stones

### UNIVERSAL / Prepolished

## 2005



L mm 9.5  
Ø mm 6.5



653.204.030.513

RA 2005

## 2006



L mm 7.5  
Ø mm 3



653.314.243.513

FG 2006



653.204.297.513

RA 2006

### UNIVERSAL / Prepolished

## 2008



L mm 10  
Ø mm 5



653.204.300.513

RA 2008

## 2012



L mm 1.5  
Ø mm 14



RA 2012






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- UNIVERSAL**
- Composite
- Ceramic
- Stones

### UNIVERSAL / Mirror effect

## 3005



L mm 9.5  
Ø mm 5.5



653.204.030.503

RA 3005

## 3006




L mm 7.5  
Ø mm 3



653.314.243.503

FG 3006



653.204.297.503

RA 3006

### UNIVERSAL / Mirror effect

## 3008



L mm 10  
Ø mm 5



653.204.300.503

RA 3008

## 3012



L mm 1.5  
Ø mm 14



RA 3012

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- COMPOSITE**
- Ceramic
- Stones

### COMPOSITE / Prepolished

## 4005



653.204.030.524

L mm 9.5  
Ø mm 6.5

RA 4005

## 4006



653.314.243.524

L mm 7.5  
Ø mm 3

FG 4006



653.204.297.524

RA 4006

### COMPOSITE / Prepolished

## 4008



653.204.300.524

L mm 10  
Ø mm 5

RA 4008

## 4012



653.204.543.524

L mm 1.5  
Ø mm 14

RA 4012

## 4014



653.204.307.524

L mm 2.5  
Ø mm 10

RA 4014

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- COMPOSITE**
- Ceramic
- Stones

### COMPOSITE / Mirror effect

## 5005



L mm 9.5  
Ø mm 6.5



653.204.030.514

RA 5005

## 5006



L mm 7.5  
Ø mm 3



653.514.243.514

FG 5006



653.204.297.514

RA 5006

### COMPOSITE / Mirror effect

## 5008



L mm 10  
Ø mm 5



653.204.300.514

RA 5008

## 5012



L mm 1.5  
Ø mm 14



653.204.543.514

RA 5012

## 5014



L mm 2.5  
Ø mm 10



653.204.804.514

RA 5014

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- COMPOSITE**
- Ceramic
- Stones

COMPOSITE / Prepolished / Diamond



### X4504



L mm 7.5  
Ø mm 6



803.204.030.523

RA X4504

### X4506



L mm 7.5  
Ø mm 3



803.314.243.523

FG X4506



803.204.297.523

RA X4506

COMPOSITE / Prepolished / Diamond



### X4508



L mm 10  
Ø mm 5



803.204.300.523

RA X4508

### X4512



L mm 1.5  
Ø mm 14



803.204.543.523

RA X4512

### X4514



L mm 2.5  
Ø mm 10



803.204.304.523

RA X4514

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- COMPOSITE**
- Ceramic
- Stones

COMPOSITE / Mirror effect / Diamond



### X5504



L mm 7.5  
Ø mm 5



803.204.030.503

RA X5504

### X5506




L mm 7.5  
Ø mm 3



803.314.243.503

FG X5506



803.204.297.503

RA X5506

COMPOSITE / Mirror effect / Diamond



### X5508



L mm 10  
Ø mm 5



803.204.300.503

RA X5508

### X5512



L mm 1.5  
Ø mm 14



803.204.543.503

RA X5512

### X5514



L mm 2.5  
Ø mm 10



803.204.304.503

RA X5514

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- Composite
- CERAMIC**
- Stones

Ceramic / Prepolished /Diamond



### X7004



L mm 7.5  
Ø mm 6



803.204.030.524

RA X7004

### X7008



L mm 10  
Ø mm 4



803.204.243.524

RA X7008

Ceramic / Prepolished /Diamond



### X7012




L mm 1.5  
Ø mm 14



803.204.543.524

RA X7012



X7008


AXIS  
dental  
SWISS QUALITY



# DENTAL BURS



Clinic



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- Composite
- CERAMIC**
- Stones

Ceramic / Mirror effect / Diamond



### X8004



803.204.030.513

L mm 7.5  
Ø mm 5

RA X8004

### X8008



803.204.243.513

L mm 10  
Ø mm 4

RA X8008

Ceramic / Mirror effect / Diamond



### X8012



803.204.543.513

L mm 1.5  
Ø mm 14


RA X8012



# DENTAL BURS



Clinic




## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- Composite
- Ceramic
- STONES

### PIERDRAS VERDE / Thin

## GR601



L mm

3

Ø 1,6 mm

19 mm

655.314.001.513

FG GR601

030

Ø 2,35 mm

22 mm

655.204.001.513

RA GR601

030

## GR638



L mm

7

Ø 1,6 mm

19 mm

655.314.110.513

FG GR638

030

Ø 2,35 mm

22 mm

655.204.110.513

RA GR638

030

### PIERDRAS VERDE / Thin

## GR645



L mm

7

Ø 1,6 mm

19 mm

655.314.161.513

FG GR645

030

Ø 2,35 mm

22 mm

655.204.161.513

RA GR645

030

## GR649



L mm

7

Ø 1,6 mm

19 mm

655.314.171.513

FG GR649

027

Ø 2,35 mm

22 mm

655.204.171.513

RA GR649

027

## GR661



L mm

6,5

Ø 1,6 mm

19 mm

655.314.286.513

FG GR661

030

Ø 2,35 mm

22 mm

655.204.286.513

RA GR661

030






# DENTAL BURS



Clinic

AXIS  
dental  
SWISS QUALITY



## DIAMOND TUGSTEN CARBIDE BURS CERAMIC POLISHERS

- Brushes
- Universal
- Composite
- Ceramic
- STONES

### Stones ARKANSAS / Extra Thin

## AR601



L mm 3

635.314.001.505

FG AR601 030



635.204.001.505

RA AR601 030

## AR638



L mm 7

635.314.110.505

FG AR638 030



635.204.110.505

RA AR638 030

### Stones ARKANSAS / Extra Thin

## AR645



L mm 7

635.314.161.505

FG AR645 030



635.204.161.505

RA AR645 030

## AR649



FG AR649 027

635.314.171.505



635.204.171.505

RA AR649 027

## AR661



L mm 6,5

635.314.268.505

FG AR661 030



635.204.268.505

RA AR661 030

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Laboratory


AXIS  
dental  
SWISS QUALITY



Our Laboratory Bur Selection offers a collection of Diamonds, Tungsten Carbide Burs and Polishers bur varieties that cover all the needs of the professional dental technician.



## How to place an order



1

### With Axis reference

Adding the diameter to the indicated Axis reference.

HP 863.016

2

### Through our e-shop

Seleccionando directamente el diámetro deseado de la punta en el e-shop. El sistema automáticamente sugerirá la referencia Axis.

HP F 863.012

3

### With ISO's number

Adding the diameter of the indicated ISO number

806.313.001.524.010






# DENTAL BURS



Laboratory



## DIAMOND

TUNGSTEN CARBIDE BURS

POLISHERS

DIVERSE

### STANDARD

Super Flex discs  
Diamond DLC  
Diverse

### Standard



### Standard

## 801



L mm



806.104.001.524

HP 801

806.104.001.534

HP G 801

## 825



L mm



806.104.304.524

HP 825

## 835



L mm



806.104.107.524

HP 835




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



## DIAMOND

### STANDARD

Super Flex discs  
Diamond DLC  
Diverse

### Standard

## 837



L mm 7



806.104.111.524

HP 837 016

## 840



L mm 9 7



806.104.111.524

HP 840 040 050

## 848



L mm 10 9



806.104.173.524

HP 848 018 040

806.104.173.534

HP G 848 018 X

## 850



L mm 10 10 9



806.104.199.524

HP 850 014 025 042

## 859



L mm 10



806.104.166.524

HP 859 018

806.104.166.514

HP F 859 018


AXIS  
dental  
SWISS QUALITY



## Laboratory



AXIS  
dental  
SWISS QUALITY



# DIAMOND

## TUNGSTEN CARBIDE BURS

## POLISHERS

## DIVERSE

## STANDARD

Super Flex discs  
Diamond DLC  
Diverse


## Standard

860



	806.104.247.524	HP 850	012	014
	806.104.247.514	HP F 850	X	014

862



006.104 249.524

HP 862 018

## Standard

863



	806.104.250.524	HP 863	012	0
	806.104.250.514	HP F 863	012	1
	806.104.250.534	HP G 863	x	0

881



	806 104 141.524	HP 881	016	02
---	-----------------	--------	-----	----




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS

POLISHERS

DIVERSE

- Standard
- SUPER FLEX DISCS**
- Diamond DLC
- Diverse

### Super Flex discs



932



Ø 2,35 mm  
44,5 mm

806 104

HP 932 450

### Super Flex discs

SF914



Ø 2,35 mm  
22 mm

L mm 0,15 0,15

Ø mm 19 22

RA SF914 190 220

Ø 2,35 mm  
44,5 mm

HP SF914 190 220

SF948



Ø 2,35 mm  
22 mm

L mm 0,15 0,15

Ø mm 19 22

RA SF948 190 220

Ø 2,35 mm  
44,5 mm

HP SF948 190 220

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



## DIAMOND

TUNGSTEN CARBIDE BURS

POLISHERS

DIVERSE

- Standard
- Disques Super Flex
- DIAMOND DLC**
- Diverse

### Diamond DLC

## DLC801



L mm 1.6



806.104.001..

HP DLC801 016

## DLC850



L mm 10



806.104.199..

HP DLC850 014

### Diamond DLC

## DLC859



L mm 10



806.104.166..

HP DLC859 018

## DLC863



L mm 10



806.104.250..

HP DLC863 016

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS




Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

Tungsten Carbide Burs are the best tool for working on soft and hard materials, depending on the professional's needs, and they reduce intervention time.



## Plaster Application



CX SG

## Plaster and Resin Applications



CX G

## Applications Alloys and Precious Metals



CX

## Applications Alloys and Precious Metals



CX F

## Applications Alloys and Precious Metals



C

## Burs 4 sides



## PLASTER APPLICATION / Super thick crossed sheets

### CX SG 251



500.104.274.223

HP CX SG 251

L mm

14

050

### CX SG 351



500.104.263.223

HP CX SG 351

L mm

14

070

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Laboratory


AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

PLASTER AND RESIN APPLICATIONS / Thick cross plates

## CX G 77



L mm 12



500.104.237.220

HP CX G 77 060

## CX G 251



L mm 14



500.104.274.220

HP CX G 251 060

## CX G 351



L mm 12



500.104.263.220

HP CX G 351 060

## CX 71



L mm 5



500.104.001.190

HP CX 71 050

## CX 77



L mm 12



500.104.237.190

HP CX 77 060

## CX 78



L mm 12



500.104.257.190

HP CX 78 060


AXIS  
dental  
SWISS QUALITY



# DENTAL BURS




Laboratory




DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

APPLICATIONS ALLOYS AND PRECIOUS METALS / Standard cross plates

## CX 79



L mm
14
15



500.104.194.190 HP CX 79 040 060

## CX 251




L mm
11,5
14




500.104.274.190 HP CX 251 040 060

## CX 257-364R




L mm
17-13




500.104.274.190 HP CX 257-364R 023

APPLICATIONS ALLOYS AND PRECIOUS METALS / Standard cross plates

## CX 261



L mm
14




500.104.194.190 HP CX 261 023

## CX 295




L mm
16




500.104.292.190 HP CX 295 023

## CX 364R



L mm
16



500.104.137.190 HP CX 364R 023




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

APPLICATIONS ALLOYS AND PRECIOUS METALS / Fine cross plates

## CX F 73



L mm 3



500.104.277.140

HP CX F 73 014

## CX F 77



L mm 5



500.104.237.140

HP CX F 77 023

## CX F 78




L mm 12



500.104.274.140

HP CX F 78 040

## CX F 79



L mm 16



500.104.190.140

HP CX F 79 060

## CX F 139




L mm 8



500.104.289.140

HP CX F 139 023

## CX F 251



L mm 14



500.104.274.140

HP CX F 251 060

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

APPLICATIONS ALLOYS AND PRECIOUS METALS / Fine cross plates


## CX F 129



L mm 8

HP CX F 129 023

## CX F 138



L mm 8

HP CX F 138 023


## C 71



L mm 2,7 3,1

HP C 71 027 031

## C 79



L mm 12

HP C 79 060

FRESA 4 LADOS

## 850.4



L mm 2,5

FG 850.4 012

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

## PRE-SYNTHEZIZED ZIRCONIA / Medium

### L0538



L mm 24  
Ø mm 10



658.104.275.523

HP L0538

## PRE-SYNTHEZIZED ZIRCONIA / Thin

### L0738



L mm 24  
Ø mm 10



658.104.275.513

HP L0738

## RESIN / Thick

### L1025



L mm 20  
Ø mm 10



658.104.300.533

HP L1025

### L1037



L mm 17  
Ø mm 6,5



658.104.275.533

HP L1037

### L1038



L mm 24  
Ø mm 10



658.104.274.533

HP L1038

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

RESIN / Medium

## L2025



L mm 24  
Ø mm 10



658.104.300.523

HP L2025

## L2037



L mm 17  
Ø mm 6.5



658.104.275.523

HP L2037

## L2038



L mm 24  
Ø mm 10




658.104.274.523

HP L2038

RESIN / Thin

## L3025



L mm 24  
Ø mm 10



658.104.300.513

HP L3025

## L3037



L mm 17  
Ø mm 6.5



658.104.275.513

HP L3037

## L3038



L mm 24  
Ø mm 10



658.104.274.513

HP L3038

AXIS  
dental  
SWISS QUALITY




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

Stones / Thick / Diamond



## LX7517



L mm 2  
Ø mm 12



805.104.042.535

HP LX7517

## LX7521



L mm 10  
Ø mm 4



805.104.173.535

HP LX7521

## LX7531



L mm 13  
Ø mm 5



805.104.107.535

HP LX7531

Stones / Medium / Diamond



## LX8517



L mm 2  
Ø mm 12



805.104.042.525

HP LX8517

## LX8521



L mm 10  
Ø mm 4



805.104.173.525

HP LX8521

## LX8531



L mm 13  
Ø mm 5



805.104.107.525

HP LX8531

AXIS  
dental  
SWISS QUALITY



# DENTAL BURS



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

Ceramic / Pre-polished / Diamond

## LX8013



L mm 2  
Ø mm 26



803.104.543.522

HP LX8013

## LX8016



L mm 2  
Ø mm 26



803.104.303.522

HP LX8016

## LX8017



L mm 2  
Ø mm 11



803.104.372.522

HP LX8017

## LX8024



L mm 13  
Ø mm 4



803.104.243.522

HP LX8024




# DENTAL BURS



Laboratory

AXIS  
dental  
SWISS QUALITY



DIAMOND TUNGSTEN CARBIDE BURS POLISHERS DIVERSE

Ceramic / Mirror effect / Diamond



## LX9013



L mm 2  
Ø mm 26



803.104.543.512

HP LX9013

## LX9016



L mm 2  
Ø mm 26



803.104.303.512

HP LX9016

## LX9017



L mm 11  
Ø mm 2



803.104.372.512

HP LX9017

## LX9024



L mm 13  
Ø mm 4



803.104.243.512

HP LX9024

AXIS  
dental  
SWISS QUALITY






# DENTAL BURS



Laboratory



DIAMANTE CARBURO DE TUNGSTENO PULIDORES DIVERSE

REINFORCED DISCS

DR C220

L mm0.3

Ø mm22

DR C220

DR M400

L mm0.5

Ø mm40

DR M400

REINFORCED DISCS

MANDRIN HP DISQUES

Ø 2.35 mm

44.5 mm

HP MANDRIN HP DISQUES

MANDRIN RA DISQUES

Ø 2.35 mm

22 mm

RA MANDRIN RA DISQUES






# DENTAL BURS



Laboratory



DIAMANTE CARBURO DE TUNGSTENO PULIDORES DIVERSE

## ADAPTATOR

4024



L mm 13  
Ø mm 2.35

RA 4024

4025



L mm 35  
Ø mm 2.35

HP 4025

4026



L mm 26  
Ø mm 2.35

RA L 4026

4032



L mm 10

RA 4032





# RADICULAR POSTS

AXIS  
dental  
SWISS QUALITY



*Effective  
solution*



AXIS  
dental  
SWISS QUALITY





# RADICULAR POSTS




Axis Dental has designed a state-of-the-art range of fiberglass radicular posts that provide behavior similar to that of a natural tooth. The heights and diameters are ideal, the restorations are simpler, stronger and safer.

## Effective, comfortable, durable and compatible

Dentists and specialists in the sector from all over the world make our radicular posts a pioneering solution in the market thanks to its multiple benefits.



4 references, 2 lengths, 2 diameters, solving 95% of the situations.





# RADICULAR POSTS



## FIBERGLASS / SET



### FIBERPOST SET

- 2 preparation strawberries
- 2 precision milling cutters
- 4 x 5 fiberglass poles
- Instructions

## FIBERGLASS / BLISTER



### FIBERPOST 1.0 - 10MM

5 fiberglass radicular posts  
Ø 1,00 - L10 mm



### FIBERPOST 1.0 - 12MM

5 fiberglass radicular posts  
Ø 1,00 - L12 mm

## FIBERGLASS / BLISTER



### FIBERPOST 0.8 - 10MM

5 fiberglass radicular posts  
Ø 0,80 - L10 mm



### FIBERPOST 0.8 - 12MM

5 fiberglass radicular posts  
Ø 0,80 - L12 mm





# RADICULAR POSTS



## BURS FOR RADICULAR POSTS



**CA XL 183.009**  
Post preparation drill 1.00



**CA XL 183.007**  
Post preparation drill 0,80

## BURS FOR RADICULAR POSTS



**CA XL 196.017**  
Post preparation drill 1.00



**CA XL 196.015**  
Post preparation drill 0,80





# TEETH WHITENING

AXIS  
dental  
SWISS QUALITY



*Effective  
solution*



AXIS  
dental  
SWISS QUALITY



# TEETH WHITENING



**With PURE teeth whitening treatments from Axis Dental you will discover a new smile in your patients with which to transmit happiness to those around them.**



## INNOVATIVE AND EFFECTIVE METHOD

PURE is a dental whitening recognized for its high prestige, backed by its technology and development.

The results, quickly visible, guarantee a whiter and brighter smile for a long time. PURE is a highly effective treatment, without risk and very easy to use. It contains an ideal pH, which provides excellent tolerance.

With PURE a difference in color is achieved in the teeth from the first days.

## EASY AND FAST TREATMENT

- Simple and safe professional treatment.
- Optimum and guaranteed long-term result.



*Quick  
and safe*





# TEETH WHITENING

AT CLINIC

AT HOME

ACCESSORIES



En el Dentista

AXIS  
dental  
SWISS QUALITY



It is the fastest treatment for whitening your teeth, achieving a whiter and brighter smile in just 1 hour. The product is applied directly to the teeth, completely protecting the soft tissues. Risk-free treatment that does not affect the structure of your teeth.



AXIS  
dental  
SWISS QUALITY

## KIT



### PURE POWER 38%

Tratamiento para 1 paciente  
38% de peróxido de hidrógeno  
- 1 jeringa de 2,5 ml  
- 1 protector gingival de 1,5 g  
- 1 PURE care de 3ml  
- Instrucción

## INDIVIDUAL SYRINGES



### PURE POWER 38% REFILL

1 jeringa de 2,5 ml  
Instrucciones



### PURE GINGIVA PROTECTOR

1 jeringa de 1,5 g  
Instrucciones



### PURE CARE

Remineralizar / Desensibilizar  
1 jeringa de 3 ml  
- 30% de Xilitol  
- 4,2% de nitrato de potasio  
- 1450 ppm de fluoruro de sodio  
Instrucciones



# TEETH WHITENING

AT CLINIC

AT HOME

ACCESSORIES



A domicilio

AXIS  
dental  
SWISS QUALITY



## Tratamiento flexible adaptado al paciente

Su dentista realizará una impresión de sus dientes, con la finalidad de elaborar sus ferulas personalizadas, facilitándole el mejor tratamiento e indicándole los pasos a seguir para optimizar el resultado.



## Three home treatment options

### PURE 10% CARBAMIDE PEROXIDE

It is the fastest treatment for whitening your teeth, achieving a whiter and brighter smile in just 1 hour. The product is applied directly to the teeth, completely protecting the soft tissues. Risk-free treatment that does not affect the structure of your teeth.

### PURE 16% CARBAMIDE PEROXIDE

It is the ideal formula for people who have a more intense tooth color and do not have sensitivity. A minimum application of 2 hours for 16 days.

### PURE PEN 16% CARBAMIDE PEROXIDE

It is a brush that provides a luminous and shiny effect on the teeth, in a simple and safe way, with just one 20-minute application. Allows between 15 and 20 applications.

## PEN



## KITS



### PURE 16% KIT

Basic treatment of 16 nights  
16% carbamide peroxide  
- 4 syringes of 3ml  
- 1 PURE care of 3ml  
- 2 thermoformable sheets  
- 1 box of splints  
- 1 color guide  
- Instructions



### PURE 10% KIT

Basic treatment of 16 nights  
10% carbamide peroxide  
- 4 syringes of 3ml  
- 1 PURE care of 3ml  
- 2 thermoformable sheets  
- 1 color guide  
- Instructions



### PURE PEN 16% KIT


Treatment of 15-20 applications  
16% carbamide peroxide  
- 1 PEN of 4ml  
- 1 LED lamp  
- 1 color guide  
- Instructions

## INDIVIDUAL SYRINGES



### PURE 16% REFILL

Additional treatment of 4 nights  
16% carbamide peroxide  
- 1 syringe of 3ml  
- Instructions



### PURE 10% REFILL

Additional treatment of 4 nights  
10% carbamide peroxide  
- 1 syringe of 3ml  
- Instructions



### PURE CARE

Remineralizing / desensitizing  
1 syringe of 3ml  
- 30% xylitol- 4.2% potassium nitrate  
- 1450 ppm sodium fluoride instructions

AXIS  
dental  
SWISS QUALITY





# TEETH WHITENING

EN CLÍNICA A DOMICILIO ACCESSORIES



### Application

The professional will make an impression of the patient's teeth in order to prepare their personalized splints, and will advise for the right choice of the most appropriate treatment, indicating the steps to follow to optimize the result.



### ACCESORIOS



We continue to innovate in order to improve our products and facilitate the task of dentists in the application of these highly demanded treatments. Consult our range of accessories on our WEBSITE.



**PURE TRAY BOX**  
Splint holder box



**PURE THERMO**  
10 thermoformable sheets  
- 12.7cm square  
- 1mm thick



**PURE LIP RETRACTOR**  
Self-cleaning



**BLOCKER**  
deposit resin  
- 1 syringe of 1.2 g  
- Instructions



### **AXIS Dental Sarl**

Ch. du Closalet 4  
1023 Crissier Switzerland  
Tel: +41 21 633 0 633  
Email: mail@axis-dental.ch

[www.axis-dental.com](http://www.axis-dental.com)

